

# **uINNOVATION - GLOBAL**

#### **Issue Highlights**

Deep Learning-Based Reconstruction of Coronary CT Angiography in Patients with Diverse Anatomical and Pathological Complexities

Dijia Wu, et al. Page 08 Clinical Evaluation of Advanced Algorithms for CCTA: A Focus on High BMI and Elevated Heart Rate Populations

John Osborne, et al. Page 24 A Vision for Magnetic Resonance Imaging to Assess Cardiotoxicity Expert Interview: Cardiac Imaging, Inc. Interview with Jens Huettges

Mary P. Watkins, et al. Page 47 Jens Huettges Page 63

#### Disclaimer

The articles contained in this magazine are provided solely by the authors, and the author(s) of each article appearing in this magazine is/are solely responsible for the content thereof as well as personal data, which is used anonymously or in compliance with applicable data privacy laws or regulations. United Imaging Healthcare makes no representation or warranties, expressly or impliedly, with respect to the accuracy, timeliness, reliability, legitimacy, applicability, fitness, originality, or completeness of the contents of this magazine. United Imaging Healthcare assumes no legal responsibility or liability for any error, omission, or illegality with respect to the material contained within.

All articles contained in this magazine only represent the opinions and views of the authors and do not implicitly or explicitly represent any official positions or policies, or medical opinions of United Imaging Healthcare or the institutions with which the authors are affiliated unless this is clearly specified. Discussions of any brand, services, or products in the magazine should not be construed as promotion or endorsement thereof.

Articles published in this magazine are intended to inspire further general scientific research, investigation, understanding, and discussion only and are not intended to and should not be relied upon as recommending or promoting a specific medical advice, method, diagnosis, or treatment by physicians for any particular individual, nor to replace the advice of a physician or other healthcare professional. Any individual wishing to apply the information in this magazine for the purposes of improving their own health should not do so without consulting with a qualified medical practitioner. All patients need to be treated in an individual manner by their personal medical advisors. The decision to utilize any information in this magazine is ultimately at the sole discretion of the reader, who assumes full responsibility for any and all consequences arising from such a decision. United Imaging Healthcare makes no representations or warranties with respect to any treatment, action, or application of medication or preparation by any person following the information offered or provided within or through the magazine.

The articles included in this magazine may contain work in progress, which represents ongoing research and development. Such technologies are not available for sale in China or the United States for clinical use and also may not be available for such sales in other countries.

Please note that the magazine is intended to be distributed within a limited scope.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you may contact us at: <u>compliance@united-imaging.com</u>

### Myocardial Perfusion Imaging by PET with Myocardial Blood Flow is Proving to be the Gatekeeper for Identifying Physiologic Severity of CAD by Guiding Treatment with Invasive Procedures and Revascularization – So Why is Adoption Limited?

Amanda Roby, MBA, PET, CNMT, RT(N)

#### 1. Introduction

Positron emission tomography (PET) myocardial perfusion imaging (MPI) with quantitative myocardial blood flow (MBF) has become the scientific standard for assessing coronary artery disease (CAD) severity, including diffuse and microvascular disease. It guides surgical intervention and ensures the safety of lifestyle-medical treatment. There are definitive scientific reports supporting the reliance on the modality, but many limitations on implementation persist. The wide acceptance of PET MPI MBF is hindered by factors such as the need for capital investment by hospitals and clinics, hardware and software limitations that scanner manufacturers have not addressed, gaps in technologist and physician education, and inconsistencies in MBF software.

The nuclear medicine stress test has been instrumental in helping cardiologists evaluate and monitor patient symptoms to identify which patients could benefit from invasive evaluation and therapy since the mid-70s. Thallium-201 (TI-201) was the first myocardial perfusion tracer used for stress, rest, and viability in a single injection. However, its low gamma emission results in lower image quality and higher radiation exposure compared to the Tc-99m or PET tracers that followed. It also requires a scan time of 4 to 24 hours. Hibernating myocardium is still identified with TI-201 today but is not commercially available in the United States. Cardiac MRI, PET, or dobutamine ECHO are the alternatives for assessing for hibernating myocardium, which have variable availability.

The Tc-99m perfusion tracers came to market in the 90s with a higher gamma peak and a shorter half-life than Tl-201allowing for the injection of a dose which is 10 times higher. This increases the target-to-background ratio improving image quality while reducing acquisition time. Unfortunately, the Tc-99m tracers lock into the mitochondria and do not redistribute like the sodium-potassium analog Tl-201 removing the ability to define hibernating myocardium. The introduction of Stress-first Tc-99m SPECT MPI improved efficiency by identifying normal patients quickly. However, without the ability to redistribute like Tl-201, patients with a stress perfusion defect require 12 to 24 hours complete the diagnosis of scar verses reversible perfusion defect. SPECT only boasts 65% sensitivity and 70% specificity. Diaphragmatic and breast attenuation defects also result in equivocal studies adding referrals to other modalities and increasing time to diagnosis, cost, radiation exposure, and risk of invasive procedures leading to unnecessary revascularization.

#### 2. Early PET Advances

The development of PET MPI with MBF occurred alongside SPECT, but PET lagged in acceptance due to limitations in computing technology, cost, and tracer accessibility. The commercial availability of Rb-82 generators in the early 2000s facilitated the expansion of cardiac PET MPI beyond academic facilities with cyclotrons, enabling the routine clinical acquisition of PET MBF. The 2000s also saw the conversion from PET with rotating rod sources to hybrid PET/CT systems. The replacement of rotating rods with CT reduced acquisition time, improved the guality of attenuation correction, and added information like coronary calcium burden, but it also brought additional struggles for PET MPI. Attenuation correction images acquired with rotating rod sources occur over minutes and account for the contraction and respiratory translation of the heart. CT for attenuation correction is captured over seconds and may lead misregistration with the PET emission. This misregistration causes artifacts in the PET data increasing the likelihood of false-positive findings.

Advances in software and hardware have addressed these issues, but PET MPI acquisition, reconstruction, and processing are still complicated. Physicians and technologists need a deeper understanding of PET physics, scanner correction structure, and coronary physiology to acquire, process, and report PET MPI MBF. 3D PET scanners are optimized to scan lower doses of 18-F for oncologic-focused exams or "hot spot" imaging, but cardiac imaging bases its principles on "cold spot" imaging. Cardiac PET utilizes "cold spot" imaging that requires accurate count recovery. Accuracy is degraded by inability or sub-optimal scatter, prompt gamma, and random corrections and deadtime and crystal saturation data loss. 3D PET scanners are more susceptible to these problems than 2D systems, especially with high-count, short half-life tracers like the most commonly used Rb-82. 3D systems, when acquired accurately, have higher sensitivity and improved image quality over 2D PET with lower radiation exposure. 2D PET systems are being phased out of production, so this discussion focuses on 3D PET. Scanner manufacturers should test the limitations of each iteration of the model they produce with the common commercial cardiac tracers like Rb-82, O-15, and N-13. Unfortunately, only F-18 and longlived isotopes like Ge-68 in phantoms are used for validation. Many factors can affect scanner limitations, such as detector material composition, electronic hardware, and computing power for corrections. Because of the range of configurations, it has been difficult for third-party MBF software providers to recommend protocol standards, delaying confidence in reporting MBF and creating confusion. Collectively, these limitations have delayed widespread adoption of cardiac PET MPI with MBF.

Over the past 15 years, cardiac PET has seen a boom in the market, with the volume of outpatient clinics tripling from 2010-2019. Increases in reimbursement in the US have aided this growth, rising from \$2250.50 to \$2750.50 just in 2022 to 2023 alone. Cardiac PET with Rb-82 allows for rest and stress to be completed in less than 30 minutes, improving the volume efficiency of the nuclear stress lab and offsetting the cost of PET operations. A generator also provides on-demand radionuclide for fast throughput of chest pain units, emergency departments, and hospitals. With few contraindications, cardiac PET can accommodate patients weighing up to 250kg with heart rate or rhythm abnormalities like atrial fibrillation, implanted devices, severe

coronary calcium burden, or impaired renal function. More importantly, the sensitivity and specificity of PET MPI are much higher than other modalities, with a sensitivity of 95%, specificity of 90% with high positive and negative predictive values comparable to cardiac catheterization, and an exam failure rate of less than 2%. However, cardiac PET is technically demanding and requires a high level of quality and physiologically informed interpretation to achieve similar precision.

#### 3. A Focus on Quality and Education

#### **3.1 Scanner Hardware and Software**

Quality cardiac PET imaging requires multiple components, including scanner design, data acquisition protocol, PET and attenuation correction alignment, and validated MBF software. Each PET scanner has unique limitations when imaging the high-count, short half-life tracers used for MBF. Optimizing protocol structure is the first step to ensure quantitative data accuracy and account for these limitations. The 3D PET list mode acquisition and correction applications are the drivers for the protocol structure and reconstruction timing. High count rate tracers and rapidly changing location dynamics of first-pass arterial activity imaging require corrections of the raw data to occur every 5 seconds or less for scatter, dead time, randoms, and prompt gamma.

The United Imaging 3D solid-state PET/CT Scanner corrects the PET data independently of the acquisition timing every 5 seconds to produce accurate quantification for any reconstruction timing protocol, whether static, dynamic, or gated. This decoupling the list mode data from the acquisition is a novel approach in list mode data correction. Other conventional PET/CT scanners apply raw data corrections based on acquisition timing structure and require careful reconstruction to ensure accurate quantification. Accurate quantification of Rb-82, for scanners with connected acquisitions and corrections, requires dynamic reconstruction of 5 seconds per frame or less for the arterial input phase (~2 minutes) and 10-30 seconds per frame for the uptake phase (~5 minutes). Moreover, it's crucial to understand the limitations of the crystals, electronics, and correction algorithms with high-count images. Each PET system requires specific testing and protocol evaluation with the tracer planned for patient use and high-count activity to assess its limitations and adapt protocols and dosing accordingly.

#### **3.2 Technologist Education**

Nuclear Medicine technologists perform flow on many body parts using a p-scope to assist in positioning with immediate visual display. Positioning the left ventricle can be challenging without this visual tool. Transitioning to PET/CT requires a greater knowledge of cross-sectional anatomy to position the heart with CT rather than a p-scope. PET/CT FOV ranges from 16-24cm in width. 16cm FOV scanners allow coverage for most hearts but can be difficult in patients with large cardiomyopathies. The preferred FOV for cardiac PET is 20-24cm. As reviewed above, using CT for cardiac attenuation correction incurs a high risk for misalignment. Co-registration of PET and CT is also not a skill taught in general nuclear medicine. With long-lived tracers, patient motion is typically corrected by reimaging. Mediation of the motion caused by heart contracting and breathing is not realistically possible. This inherent cardiac and respiratory motion, combined with the 75-second half-life of Rb-82, requires technologists to coach their patients to remain still during acquisition phases. Manual or automated co-registration software is available across all manufactures of modern PET/CT scanners. Registration should allow for translation in the x, y, and z-axis. United Imaging takes this further by adding a rotational shift to its alignment software [1, 2]. Every PET cardiac image must be assessed and corrected for misalignment with the attenuation map, even on systems with automated registration techniques. Any relative perfusion defect not supported by patient history, symptoms, or other findings requires technologists and physicians to scrutinize the images a second time for misalignment. Alignment of PET and attenuation maps for cardiac PET is a specialized technologist skill, not comparable to any other exam processing in nuclear medicine. Hence technologist experience and training play a large role in image quality. Cardiac PET images typically require at least two iterations of reconstruction. The first is the initial unaligned images, then repeating all reconstructions with an aligned attenuation correction. These steps are not automated within the acquisition protocols provided by any manufacture and require additional expertise and time compared to oncologic PET/CT.

Some technologist processing skills are transferable from nuclear medicine cardiac imaging, including angulation and apex, base, and contouring limits to exclude extracardiac activity and ejection fraction calculation. Performing these tasks accurately and reliably is paramount for comparing artery-specific perfusion changes from rest to stress and day to day for comparisons. Processing MBF comprises new theory, processing skills, and critical thinking outside the skill set of an entry level technologist. Many MBF software platforms are available for purchase but vary greatly in concept and accuracy. Johnson's 2021 editorial summarizes how buyers should compare each software vendor to a comprehensive list of validation criteria [3]. The list includes validation in animal models, normal volunteers to ischemic and infarcted patients, precision established by test-retest, customizable arterial input, and documented for clinical utility by use in a clinical outcome review, economic benefit, and randomized trial.

Once an institution decides on a flow package, technologists must be trained to process the data for visualization, flow, and assessment for quality. Unfortunately, current nuclear medicine training programs have neither the clinical sites nor time required to train students in cardiac PET. Cardiac PET is a specialty rotation, allowing a maximum of 4 to 5 weeks of instruction. Our institution's experience in cross-training nuclear medicine technologists in cardiac PET suggests a minimum training period of 3 to 6 months to gain competency. The time to train a technologist depends on the complexity of the PET/CT system they are using and their base knowledge of Rb-82 imaging. To address the workforce education gap, the University of Texas Medical School -Houston has created a Cardiac PET Workforce Training program. Post-graduate nuclear medicine technologists spend six months gaining a comprehensive education in cardiac PET/CT exams, including stress testing, viability, and inflammation imaging as well as certification on two different brands of PET/CT systems, the Cardiogen Rb-82 infusion system, and HeartSee quantitative perfusion software.

Adding flow acquisition and processing to a facility that is currently acquiring Rb-82 presents a smaller learning curve than creating a new service line of cardiac PET with MBF. Cardiac software packages attempt many tasks automatically for the technologist, which creates efficiency for experienced technologists but can degrade quality for inexperienced technologists. Without proper training, technologists will lack the ability to connect patient history, stress results, and image findings to assess accuracy.

Time-activity curve (TAC) modeled software has automated alignment that assists with alignment of the PET flow phase

(24+ frames) and uptake (10+ frames) within segmental cardiac boundaries but is not routinely complete. Misalignment causes substantial errors in quantitative perfusion and technologists must check, correct when necessary, and recognize when those boundaries, though aligned, cannot accurately define the TAC modeled MBF due to the rigid segmentation [4].These tasks are typically considered physician-level duties, but unfortunately, physicians are also not trained to process PET MBF in radiology or cardiology fellowship programs to address this issue, the UT Cardiology fellowship program has recently added a one-month rotation in cardiac PET.

Frame alignment is unnecessary for compartmental model packages that do not use segmental boundaries for the flow and myocardium. Removing segmental boundaries also improves accuracy as averaged segmental flow values in a 17segment model that arbitrarily overlap arterial distributions are not reported. The overlap causes confusion, routinely under-representing severity, when comparing PET MBF to IAC-FFR. PET's high resolution allows for pixel-level expression of severity, substantially increasing accuracy of an artery-specific diagnosis over the aged 17-segment model carried over from nuclear medicine SPECT [5]. Furthermore, third-party cardiac software often only allows for visualization of the final MBF processing, not the console reconstruction and alignment. Additionally, Radiology Business magazine recently reported that "41% of radiologists surveyed expressed doubt in their skills processing cardiac images," along with 28% doubting their ability in nuclear medicine, despite documented competence [6]. The cardiac PET technologist plays an integral role by applying their knowledge of nuclear physics, instrumentation, coronary physiology, and each patient's unique history to present the most accurate data to the physician for interpretation while working in high-throughput clinics.

#### **3.3 Clinical Gatekeeper**

The Supporting Science: Over extensive literature publication for the past 40 years, the University of Texas cardiac PET Center has established the physiologic basis, technological requirements, and clinical application of PET using HeartSee perfusion software. Coronary Flow Capacity (CFC) maps by HeartSee offer particular value to invasive cardiology not provided by any other invasive or non-invasive technology [7]. Coronary angiographic anatomy measurement of percent stenosis fails to adequately stratify CAD risk severity because risk is derived from coronary blood flow and it depends on the arterial lumen radius raised to the 4th power. Therefore, the angiogram's limited resolution cannot accurately measure small changes in stenosis dimension that may cause life threatening changes. The arteriogram is commonly considered the "gold standard" of CAD, but extensive literature proves otherwise. Every randomized trial of coronary stents or bypass surgery guided by angiographic severity has failed to reduce mortality.

Invasive simulated fractional flow reserve (IAC-FFR) based on intra-coronary pressure during maximal coronary blood flow has advanced beyond the angiogram for assessing coronary stenosis. However, it remains inferior to CFC by PET for two reasons. Randomized trials of FFR guided revascularization have failed to demonstrate a reduction in mortality since IAC-FFR also fails to account for coronary blood flow or myocardial perfusion. Cardiac PET does so accurately and quantifies mortality risk and its predicted improvement after revascularization.

PET CFC identifies patients with high mortality that reduces by 54% after PET guided revascularization compared to norevascularization for comparable PET severity in large nonrandomized cohorts of 7000 patients followed for 14 years for major adverse coronary events (MACE) [8-13]. Figure 1 summarizes the transformation of a single observation in a single artery to not only quantifying perfusion per pixel, but the impact of revascularization related to flow by PET taken from the Gould et al publication in 2022. The PET guided survival benefit has now been confirmed by the randomized CENTURY trial being reported at the European Society of Cardiology (ESC) congress in August 2023. The recently completed randomized CENTURY trial also indicates that clinical management and interventions guided by PET HeartSee perfusion software reduced death, myocardial infarction, late revascularization and major adverse coronary events (MACE) compared to standard community care without PET. The HeartSee CFC display condenses the complexity of PET perfusion and MBF giving physicians and patients reassurance to safely pursue lifestyle-medical management and reserving invasive interventions only for objectively severe, high-risk, perfusion abnormalities.

PACIFIC and ReASSESS trials show that PET MBF has the highest accuracy rate (86%) compared to intraarterial cardiac catheterization fractional flow reserve (IAC-FFR), verses SPECT or FFRCT accuracy rates of 68-76% and 70%,

respectively. PET also accommodates patients with atrial fibrillation, metal implants, high BMI, and dense coronary calcium, unlike FFRCT and cardiac MRI (cMRI). And although technically difficult, PET MPI with MBF reports a much lower exam failure rate of 2% compared to FFRCT of 13%[14-16] If the goal is patient-centered care, we must choose the exam that provides the patient with the most accurate and comprehensive result.

A recent case represents the limited scope of anatomical imaging. A 66-year-old female presented to the ER with chest and arm pain radiating to shoulder and back accompanied with vomiting. The day prior she was found to have COVID after feeling generalized malaise and unwell for a week. She had a past history of CAD and PCI to LAD in 2019, anomalous origin of the RCA from the left coronary sinus with course between the aorta and pulmonary artery, as well as hypertension, hyperlipidemia, diabetes, and obesity. Her high-sensitivity troponins peaked at ~1800 and downtrended, anterolateral T-wave inversion without ST-segment elevation. The outlying hospital diagnosed NSTEMI and transferred her for possible intervention of the anomalous RCA.

The patient underwent a left heart catheterization (LHC) which showed a patent stent to the LAD, no stenosis of the left circumflex, but unable to engage the RCA. Cardiac CTA confirmed anomalous origin of a right dominant system from the left coronary cusp, with an intramural course and mild to moderate stenosis with calcified atherosclerotic plaque. Two days later, the patient underwent a second LHC which successfully engaged a patent RCA. The care team considered RCA unroofing verses single vessel CABG but ordered quantitative rest-stress PET perfusion imaging to confirm RCA as source of angina. Dipyridamole stress replicated her angina and caused significant ST depression despite PET MPI showing no resting or stress induced perfusion defect and good CFC with normal regional and global rest-stress flows of 1.76 and 2.23 cc/min/g respectively.

Relative tomographic images revealed reduced subendocardial to subepicardial perfusion ratio due to left ventricular hypertrophy and diffuse non-obstructive epicardial coronary atherosclerosis during high coronary flow.

Figure 2 displays this patient's course of events that finally achieved an answer for her angina and a plan for her care. The PET report recommended vigorous blood pressure and heart rate control (201/80 mmHg and 81 bpm at baseline) to decrease the myocardial demand causing reduced hyperemic coronary pressure and subendocardial ischemia. Slowing the heart rate also provides longer diastolic perfusion time for better perfusion through a thickened LV wall. Prior to PET, the invasive and non-invasive anatomical imaging failed to identify the cause of this patient's chest pain after 9 days of in-patient care, accumulating 2215 mGy of radiation exposure and administering 330 ml of contrast media. The PET MPI with MBF provided the care team and the patient confidence to cancel the planned surgery on the anomalous RCA in favor of medical management.

#### 3.4 Summary

While substantial capital investment is required to initiate a PET MPI with MBF program, its downstream efficiency compensates by reducing expenses throughout the care system from less accurate modalities and their associated costs and lost time. It is accepted that accurately quantified cardiac PET is the most robust tool for assessing the severity and directing clinical management of the CAD patient. Adopting such imaging protocol confidently requires investing personnel time and resources to understand the PET/CT system, quantitative perfusion software, coronary pathophysiology and its clinical application. Achieving the best patient care at an efficient cost necessitates collaboration between physicians and technologists unmatched by any other invasive or non-invasive technology.



#### From Concept to Clinical Coronary Pathophysiology

Figure 1 Quantifying perfusion per pixel, and the impact of revascularization related to flow by PET



Figure 2: Case study of 66-year-old female presenting with chest pain and complex coronary angiogram for which planned bypass surgery was avoided by quantitative reststress PET perfusion imaging (see text for details). (A) LHC patent LAD, no stenosis of LCx, unable to engage the RCA. (B) Cardiac CTA suggesting moderate narrowing of the RCA. (C) repeat LHC found RCA to have mild to moderate stenosis. (D) Reduced stress relative subendo/subepicardial ratio by HeartSee software. (E) HeartSee PET rest-stress relative perfusion and CFC maps. (F) Tomographic PET images. (G) ECG at admission with T-wave inversion.

#### References

- Loghin C, Sdringola S, Gould KL. Common artifacts in PET myocardial perfusion images due to attenuation-emission misregistration: clinical significance, causes, and solutions. J Nucl Med. 2004 Jun;45(6):1029-39. PMID: 15181138.
- 2. Slomka P, Diaz-Zamudio M, Dey D, Motwani M, Brodov Y, Choi D, Hayes S, Thomson L, Friedman J, Germano G, Berman D. Automatic registration of misaligned CT attenuation correction maps in Rb-82 PET/CT improves detection of angiographically significant coronary artery disease. doi:10.1007/s12350-014-0060-9
- 3.Johnson NP, Gould KL. How shall we judge a PET flow model? J Nucl Cardiol. 2022 Oct;29(5):2551-2554. doi: 10.1007/s12350-021-02805-5. Epub 2021 Sep 24. PMID: 34561847
- 4. Koenders SS, van Dijk JD, Jager PL, Ottervanger JP, Slump CH, van Dalen JA. How to detect and correct myocardial creep in myocardial perfusion imaging using Rubidium-82 PET? J Nucl Cardiol. 2019 Jun;26(3):729-734. doi: 10.1007/s12350-019-01650-x. Epub 2019 Feb 20. PMID: 30788759; PMCID: PMC6517341.
- 5. Johnson NP, Gould KL Retention models: 'tis the gift to be simple. J Nucl Cardiol 2022; 5:2595-2598.
- 6. Radiologists read across an average of 5 subspecialties but aren't always confident doing so. June 22, 2023. https://radiologybusiness.com/topics/healthcaremanagement/medical-practice-management/radiologistsread-across-average-5-subspecialties-arent-alwaysconfident-doing-so. Accessed July 20, 2023.
- 7. Johnson NP, Gould KL. Physiologic basis for angina and ST change: PET-verified thresholds of quantitative stress myocardial perfusion and coronary flow reserve. J Am Coll Cardiol: Cardiovascular Imaging 2011;4:990-998. Awarded the 2011 Young Author Achievement Award by the Journal of the American College of Cardiology Cardiovascular Imaging.
- B. Gould KL, Nguyen T, Kirkeeide RL, Roby AE, Bui L, Kitkungvan D, Patel MB, Madjid M, Haynie MP, Lai D, Li R, Narula J. Subendocardial and Transmural Myocardial Ischemia: Clinical Characteristics, Prevalence and

Outcomes With and Without Revascularization. J Am Col Cardiol CV Imaging 2023;16:78-94.

- 9. Gould KL, Kitkungvan D, Johnson NP, et al. Mortality Prediction By Quantitative PET Perfusion Expressed As Coronary Flow Capacity With and Without Revascularization. J Am Coll Cardiol CV Imaging 2021;14:1020-34.
- 10.Gould KL, Nguyen TT, Kirkeeide R, Johnson NP. Atlas of Nuclear Cardiology, Chapter 6 Clinical Coronary Physiology and Quantitative Myocardial Perfusion, Fifth Edition, Editors Dilsizian and Narula, 2021, Springer Nature Switzerland.
- 11.Gould KL, Johnson NP, Roby AE, et al. Regional artery specific thresholds of quantitative myocardial perfusion by pet associated with reduced mi and death after revascularization in stable coronary artery disease. J Nucl Med 2019;60:410-417.
- 12.Gould KL, Johnson NP, Bateman TM, et al. Anatomic versus Physiologic Assessment of Coronary Artery Disease: Role of CFR, FFR, and PET Imaging in Revascularization Decision-Making. J Am Coll Cardiol 2013;62:1639–53.
- 13.Gould K, Johnson NP. Coronary physiology: Beyond CFR in microvascular angina. J Am Coll Cardiol 2018;72:2642-62 (Selected for Highlights of the year 2018 in J Am Coll Cardiol 2019;73:849).
- 14.Driessen RS, Danad I, Stuijfzand WJ, Raijmakers PG, Schumacher SP, van Diemen PA, Leipsic JA, Knuuti J, Underwood SR, van de Ven PM, van Rossum AC, Taylor CA, Knaapen P. Comparison of Coronary Computed Tomography Angiography, Fractional Flow Reserve, and Perfusion Imaging for Ischemia Diagnosis. J Am Coll Cardiol. 2019 Jan 22;73(2):161-173. doi: 10.1016/j.jacc.2018.10.056. PMID: 30654888.
- 15.Sand NPR, Veien KT, Nielsen SS, Nørgaard BL, Larsen P, Johansen A, Hess S, Deibjerg L, Husain M, Junker A, Thomsen KK, Rohold A, Jensen LO. Prospective Comparison of FFR Derived From Coronary CT Angiography With SPECT Perfusion Imaging in Stable Coronary Artery Disease: The ReASSESS Study. JACC Cardiovasc Imaging. 2018 Nov;11(11):1640-1650. doi: 10.1016/j.jcmg.2018.05.004. Epub 2018 Jun 13. PMID: 29909103.

16.Johnson NP, Gould KL. Coronary physiology – Simulations can't beat the real thing. J Am Coll Cardiol CV Imaging 2020 (In press on line, https://doi.org/10.1016/j.jcmg.2020.02.014)

#### **Author Biography**



Amanda E Roby, PET, CMNT, RT(N), MBA Weatherhead PET Imaging Center Assistant Director PET Technical Operations

Amanda Roby is a graduate of the University of Oklahoma Nuclear Medicine program with 18 years of experience predominately in cardiac PET, and Assistant Director of Technical Operations at the University of Texas Weatherhead PET/CT Imaging Center. The Weatherhead PET Center is located in the heart of the Texas Medical Center. She worked with the team to invent and validate HeartSee, a myocardial perfusion and quantitative coronary flow software, and the measurement of Coronary Flow Capacity. She has contributed many publications on the subject to their extensive databas of over 12000 PETs and 15 years of patient follow up data.

## **PASSION for CHANGE**

©2022 United Imaging Healthcare Co., Ltd. All rights reserved.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you are welcome to contact us at the following email address: <u>uinnovation-global@united-imaging.com</u>