

u!nnovation.

**UNITED
IMAGING** 



Global Spread of external contributors



Foreword

A Note from the Editorial Board

Welcome to this year's annual magazine from the Global Scientific Collaboration team.

uInnovation is a scientific magazine published by United Imaging Healthcare that has been successfully distributed for over past three years. It aims to serve as a platform for sharing ground-breaking advancements, emerging trends, and future possibilities in the vast expanse that is oncology.

uInnovation is currently in its fourth edition. This year's edition will inform, engage, and inspire you about the latest developments and applications of United Imaging Healthcare. This journal includes quick read sections for those in a rush, and appealing images to promote visual understanding.

United Imaging Healthcare is a provider of high-end medical equipment and medical IT solutions. From our headquarters in Shanghai's Jiading district to our network of research and development centers throughout the world, our global mission is to provide medical institutions with a full range of healthcare solutions, from diagnostic imaging and radiation therapy equipment to service, training, and medical IT solutions.

We are dedicated to expanding access to high-quality medical care and improving the value of our services.

VISION — Leading Healthcare Innovation

MISSION — To Bring Equal Healthcare for All

To learn more, visit <https://www.united-imaging.com/>

Disclaimer

The articles contained in this magazine are provided solely by the authors, and the author(s) of each article appearing in this magazine is/are solely responsible for the content thereof as well as personal data, which is used anonymously or complied with applicable data privacy laws or regulations. United Imaging Healthcare makes no representation or warranties, expressly or impliedly, with respect to the accuracy, timeliness, reliability, legitimacy, applicability, fitness, originality, or completeness of the contents of this magazine. United Imaging Healthcare assumes no legal responsibility or liability for any error, omission, or illegality with respect to the material contained within.

All articles contained in this magazine only represent the opinions and views of the authors and do not implicitly or explicitly represent any official positions or policies, or medical opinions of United Imaging Healthcare or the institutions with which the authors are affiliated unless this is clearly specified. Discussions of any brand, services, or products in the magazine should not be construed as promotion or endorsement thereof.

Articles published in this magazine are intended to inspire further general scientific research, investigation, understanding, and discussion only and are NOT intended to and should not be relied upon as recommending or promoting a specific medical advice, method, diagnosis, or treatment by physicians for any particular individual, nor to replace the advice of a medical doctor or other healthcare professional. Any individual wishing to apply the information in this magazine for the purposes of improving their own health should not do so without consulting with a qualified medical practitioner. All patients need to be treated in an individual manner by their personal medical advisors. The decision to utilize any information in this magazine is ultimately at the sole discretion of the reader, who assumes full responsibility for any and all consequences arising from such a decision. United Imaging Healthcare makes no representations or warranties with respect to any treatment, action, or application of medication or preparation by any person following the information offered or provided within or through the magazine. United Imaging Healthcare shall remain free of any fault, liability, or responsibility for any loss or harm, whether real or perceived, resulting from the use of information in this magazine.

The articles included in this magazine may contain work in progress, which represents ongoing research and development. Such technologies are not available for sale in the United States for clinical use and also may not be available for such sales in other countries around the world.

Please note that the magazine is intended to be distributed only within a limited scope instead of publication.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you are welcomed to contact us at the following email address: compliance@united-imaging.com

Total-Body PET/CT in Theranostics: Advancing Molecular Imaging in Precision Oncology

Harshad R. Kulkarni

BAMF Health, Grand Rapids, Michigan, UNITED STATES

Abstract

Total-body PET/CT is emerging as a transformative technology in theranostics—the integration of diagnostic imaging with targeted radionuclide therapy. By employing long axial fields-of-view (FOV) PET/CT and increased sensitivity, these scanners are reshaping oncology imaging and patient management. This article introduces theranostics and total-body PET/CT, reviews recent advances, and discusses their clinical impact through real-world experience. The unprecedented sensitivity and speed of total-body imaging enhance lesion detection, therapy planning, and treatment monitoring, while enabling more precise and personalized cancer care. Key challenges, including cost and data management, are addressed alongside future directions in dosimetry and novel radiopharmaceutical development.

" Whole-body PET/CT technology isn't just about pixels or counts; it's about giving patients a fighting chance where cancer has nowhere to hide."

Theranostics

Over the past two decades, this “see what you treat” strategy has been firmly established as a pillar of precision oncology (1). Foundational clinical work in neuroendocrine tumors (NETs) demonstrated the feasibility and efficacy of peptide receptor radionuclide therapy (PRRT), and early experience with prostate-specific membrane antigen (PSMA)-directed radioligand therapy broadened the scope of theranostics in prostate cancer (2,3). The advent of PSMA PET imaging provided the crucial diagnostic backbone, enabling reliable target visualization, patient selection, and therapy monitoring (4). Together, these foundational efforts laid the basis for routine adoption of PET-guided therapy selection and monitoring in oncology (5,6).

Building on this foundation, landmark phase 3 trials have confirmed the clinical value of theranostics. In NET, the NETTER-1 trial showed that [177Lu] Lu-DOTA-TATE (177Lu-

Dotatate) significantly improved progression-free survival (PFS), leading to regulatory approval, and the recent NETTER-2 trial demonstrated benefit in higher-grade gastroenteropancreatic NETs (7,8). In prostate cancer, the pivotal VISION trial showed that [177Lu]Lu-PSMA-617 (177Lu-PSMA-617) improved overall survival (OS) in patients with advanced castration-resistant disease, while the PSMAfore trial moved PSMA radioligand therapy earlier into the taxane-naïve metastatic castration-resistant prostate cancer (mCRPC) setting with radiographic PFS benefit (9,10). Looking ahead, the multicenter WARMTH Act study underscored the promise of 225Ac-PSMA as an alpha-emitter option after progressing on 177Lu-PSMA (11). Across these studies, PET serves as a gatekeeper for target expression and as a tool for serial response assessment, cementing its role in theranostics workflows.

Total-Body PET/CT

The concept of an extended axial FOV PET arose to the need to maximize sensitivity and enabling comprehensive kinetic imaging across distributed diseases. In 2018, Cherry and colleagues reported the first prototype and early human studies, demonstrating that with total-body PET/CT entire body could be scanned in one bed position and it could improve the detection of avid tracer uptake while shortening acquisition time (12). Traditional PET/CT systems cover ~20–30 cm per bed, requiring multiple passes for whole-body coverage. In contrast, total-body PET/CT systems extend to ~1–2 meters, enabling true head-to-toe imaging in a single acquisition. Depending on AFOV length, detector technology, and reconstruction, effective sensitivity gains on the order of ~10–40× have been reported, permitting high-quality scans in minutes (12).

One immediate benefit is enhanced lesion detectability. Long axial FOV PET improves the detection of small lesions and low-contrast lesions compared with conventional scanners (13). The ability to lower administered activity while maintaining non-inferior diagnostic performance has been validated in selected settings, supporting dose reduction and

patient-centric protocols (14). Furthermore, advanced imaging protocols on total-body PET/CT—dynamic acquisition and delayed imaging –improving contrast-to-noise and quantitative accuracy beyond static imaging alone (15). These advantages are now being demonstrated in prospective clinical settings. Early studies have shown that total-body PET can achieve diagnostic-quality images in a fraction of the time traditionally required, underscoring the potential for dramatically faster workflows and greater

patient comfort (16). Simultaneously, total-body PET/CT systems facilitate first-in-human evaluations of novel radiopharmaceuticals by tracking whole-body distribution and pharmacokinetics in unprecedented detail (17). Together, these applications refine established protocols that illustrate and accelerate the development and translation of new theranostic agents—bridging discovery and clinical care with speed, sensitivity, and scalability.




Fig. 1 A 65-year-old man with Gleason 4+5 metastatic prostate cancer previously treated with nine cycles of ¹⁷⁷Lu-PSMA radioligand therapy at an outside center. He initially demonstrated a good response, followed by continued disease progression documented on serial ⁶⁸Ga-PSMA-11 total-body PET/CT (A–D, sequential MIP images).

In theranostics, sensitive detection, serial monitoring, and pharmacokinetics /dosimetry are crucial—these advances are transformative. Total-body PET/CT enhances patient care through earlier detection and individualized monitoring, and also opens pathways for innovation in tracer development, dosimetry, and adaptive therapy strategies.

Scaling Theranostics with Total-Body PET/CT:
Impact, Challenges, and Future Potential
Speed, sensitivity, and quantitation

Total-body PET/CT has reshaped theranostic practice by bringing speed, sensitivity, and quantitation into daily practice. Exams that once required 30–40 minutes can often be completed within 5 minutes, mitigating patient burden, motion artifacts, and scheduling bottlenecks, thereby enabling a significantly higher patient throughput. Higher sensitivity captures subtle disease—tiny nodes, early marrow/bone lesions, or faint visceral deposits—that directly influence eligibility and treatment planning. Total-body dynamic acquisitions add another dimension: pharmacokinetics and dosimetry may be measured in a

single visit, shifting the paradigm from “image and infer” to “measure and personalize”. In routine operations, where full research dosimetry is not yet standard or reimbursed,

protocolized quantitative PET and safety labs care are increasingly evaluated in trials and select clinical pathways.




Fig. 2 An 81-year-old man with stage IVB (cTXN0M1b, Gleason 4+4=8, Grade Group 4) prostate adenocarcinoma (pre-treatment PSA 56 ng/mL). He received androgen deprivation therapy followed by apalutamide, pelvic external beam radiation, and 5 cycles of docetaxel chemotherapy (nadir PSA 1.9 ng/mL), complicated by congestive heart failure and peripheral neuropathy. (A-D) Pre-therapy ^{68}Ga -PSMA-11 PET/CT demonstrates extensive PSMA-avid osseous metastases, predominantly in the left femur, tibia, and tarsal bones, as well as additional lesions in the right lower limb and pelvis.

A scalable care model

At BAMF Health, these capabilities are embedded in an end-to-end pathway: standardized referral and triage, on-site GMP radio-pharmacy support, target-specific PET for eligibility, multidisciplinary review, therapy delivery, and post-therapy molecular response assessment. High-throughput imaging supports high-throughput, data-informed care. Robust molecular responders may be de-escalated or paused to limit exposure, and treatment continuation is guided by standardized imaging and clinical criteria. Throughout, standardized acquisition, reconstruction, and calibration and quality assurance (QA) maintain quantitative rigor as volumes grow.

PET/CT plus PET/MR

Co-location with PET/MR adds complementary strengths. PET/CT offers speed, sensitivity, and throughput; PET/MR provides superior soft-tissue contrast, diffusion metrics, and motion-robust sequences. Hybrid strategies include rapid total-body PET/CT for eligibility and quantity followed by targeted PET/MR for problem regions, and PET/MR for organ-specific monitoring when minimizing radiation is desirable. Together, the modalities support precise target confirmation, confident lesion detection, and cleaner longitudinal response curves.




Fig. 3 A patient with Gleason 4+3=7 prostate adenocarcinoma after androgen deprivation therapy, novel androgen-axis therapy, and docetaxel. ^{68}Ga -PSMA-11 PET/CT demonstrates widespread PSMA-positive lymph node and osseous metastases. The patient was deemed eligible for ^{177}Lu -PSMA radioligand therapy. (A-D) MIP images with varying acquisition times: (A) 10 seconds, (B) 30 seconds, (C) 2 minutes, (D) 5 minutes. Even the 30-second acquisition clearly delineates the PSMA-positive disease burden.

Refining therapy response

Total-body PET/CT enables more consistent quantification of whole-body tumor burden by combining complete coverage with high sensitivity. With complete body coverage

and superior sensitivity, cumulative uptake metrics can be quantified more reliably, enabling consistent tracking of global tumor burden. In NETs, reductions in somatostatin receptor-positive tumor volume (SSTR-TV) have been shown

to align with clinical improvement after PRRT. In prostate cancer, accurate quantification on PSMA PET enables whole-body tumor burden measurement and offers prognostic insights that can help guide decision-making. Beyond static measures, dynamic and longitudinal imaging provides access to biomarkers that may predict treatment efficacy early in the course of therapy. Faster, lower-dose acquisitions also open the door to more frequent monitoring, enabling nimble, evidence-driven adjustments.

Remaining challenges

Capital expenditure and facility build-out limit the access of LAFOV systems. Dynamic acquisitions generate vast data volumes, demanding robust PACS/VNA infrastructure and

GPU-enabled analytics. Multicenter harmonization requires phantom experiment, cross-calibration, and locked or standardized reconstructions to preserve SUV fidelity. Workforce readiness—from physicians and physicists to technologists and nursing—must keep pace via standardized training. Supply-chain fragility for isotopes and reimbursement gaps for dosimetry and dynamic studies also constrain scalability. Finally, interpretation requires discipline, as physiological uptake variants are more conspicuous on highly sensitive images.




Fig. 4 ⁶⁸Ga-DOTATATE total-body PET/CT (A, MIP) in a patient with resected, well-differentiated ileal neuroendocrine tumor (NET) revealed a solitary liver lesion in segment VIII. Subsequent PET/MR performed with the same injection confirmed the lesion as a somatostatin receptor-positive metastasis (B, MRI; C, fused PET/MR).

The near future

Looking ahead, several trajectories are within reach. Routine prospective dosimetry with short multi-time-point scans is likely to become standard for PRRT and PSMA radiopharmaceutical therapies, improving dose-response correlations. Pre-therapeutic dosimetry using longer-lived

positron emitters (e.g., ⁶⁴Cu- or ⁸⁹Zr-labeled analogues/antibodies) allows imaging over extended periods to model tumor and organ kinetics before therapy. This provides a pathway toward true individualized therapy planning in precision oncology. AI-assisted quantitation could reduce reporting time and

inter-reader variability by automating lesion detection, whole-body tumor burden metrics, and kinetic mapping. Fast acquisitions further support the concept of a “one-stop” therapy day: eligibility PET in the morning, therapy midday, and a short post-therapy verification scan in the afternoon. Expanding tracer portfolios, including FAP, GRPR, CXCR4, HER2, and antibody-based agents—can be evaluated rapidly in early-phase studies, accelerating drug development and translation. Finally, networked theranostics built on harmonized imaging standards will allow patients to be imaged and treated closer to home. In this vision, high-sensitivity long axial FOV PET/CT systems serve as the enabling technology, linking advanced quantitation with scalable, patient-centered care.

Why does this matter?

For patients, these advances translate into shorter and less stressful scans, earlier detection of actionable disease, and faster therapy adjustments when treatments underperform. For health systems, throughput and quantitation expand capacity while raising the bar for precision. And for innovators, total-body PET/CT provides the measurement engine that turns theranostics from promise into reproducible, scalable care—the very model BAMF is designed to deliver.

References

1. Baum RP, Kulkarni HR. (2012) THERANOSTICS: From Molecular Imaging Using Ga-68 Labeled Tracers and PET/CT to Personalized Radionuclide Therapy - The Bad Berka Experience. *Theranostics*.
2. Kwekkeboom DJ, et al (2008) Treatment with the radiolabeled somatostatin analog [177 Lu-DOTA 0,Tyr3]octreotate: toxicity, efficacy, and survival. *J Clin Oncol*.
3. Kulkarni HR, et al (2016). PSMA-Based Radioligand Therapy for Metastatic Castration-Resistant Prostate Cancer: The Bad Berka Experience Since 2013. *J Nucl Med*.
4. Afshar-Oromieh A, et al (2013). PET imaging with a [68Ga] gallium-labelled PSMA ligand for the diagnosis of prostate cancer: biodistribution in humans and first evaluation of tumour lesions. *Eur J Nucl Med Mol Imaging*.
5. Kulkarni HR, Baum RP. Patient selection for personalized peptide receptor radionuclide therapy using Ga-68 somatostatin receptor PET/CT. *PET Clin*. 2014
6. Kulkarni HR, Baum RP. Theranostics with Ga-68 somatostatin receptor PET/CT: monitoring response to peptide receptor radionuclide therapy. *PET Clin*. 2014 Jan;9(1):91-7.
7. Strosberg J, El-Haddad G, et al (2017) NETTER-1 Trial Investigators. Phase 3 Trial of 177Lu-Dotatate for Midgut Neuroendocrine Tumors. *N Engl J Med*.
8. Singh S, et al (2024); all the NETTER-2 Trial Investigators. [177Lu]Lu-DOTA-TATE plus long-acting octreotide versus high dose long-acting octreotide for the treatment of newly diagnosed, advanced grade 2-3, well-differentiated, gastroenteropancreatic neuroendocrine tumours (NETTER-2): an open-label, randomised, phase 3 study. *Lancet*.
9. Sartor O, et al (2021); VISION Investigators. Lutetium-177-PSMA-617 for Metastatic Castration-Resistant Prostate Cancer. *N Engl J Med*.
10. Fizazi K, Chi KN, et al (2025); PSMAfore Investigators. Final overall survival and safety analyses of the phase III PSMAfore trial of [177Lu]Lu-PSMA-617 versus change of androgen receptor pathway inhibitor in taxane-naïve patients with metastatic castration-resistant prostate cancer. *Ann Oncol*.
11. Sathekge MM, et al (2024) Actinium-225-PSMA radioligand therapy of metastatic castration-resistant prostate cancer (WARMTH Act): a multicentre, retrospective study. *Lancet Oncol*.
12. Cherry SR, et al (2018). Total-Body PET: Maximizing Sensitivity to Create New Opportunities for Clinical Research and Patient Care. *J Nucl Med*.
13. Viswanath V, et al (2020) Numerical observer study of lesion detectability for a long axial field-of-view whole-body PET imager using the PennPET Explorer. *Phys Med Biol*.
14. Sachpekidis C, et al (2023) Application of the long axial field-of-view PET/CT with low-dose [18F]FDG in melanoma. *Eur J Nucl Med Mol Imaging*.
15. Wu Y, et al (2024) The role of dynamic, static, and delayed total-body PET imaging in the detection and differential diagnosis of oncological lesions. *Cancer Imaging*.
16. Kulkarni HR, et al (2023) *J Nucl Med* 2023;64 (supplement 1): P1408 https://jnm.snmjournals.org/content/64/supplement_1/P1408
17. Kulkarni HR, et al (2024) First-in-Human Total-Body PET/CT Imaging Using 89Zr-Labeled MUC5AC Antibody in a Patient with Pancreatic Adenocarcinoma. *J Nucl Med*.

Author Biography



Dr Harshad R. Kulkarni, MD
Chief Medical Advisor, BAMF Health,
Grand Rapids, MI, UNITED STATES

Harshad R. Kulkarni, MD, is an internationally recognized physician-scientist and Chief Medical Advisor at BAMF Health (Grand Rapids, MI), where he oversees strategy for clinical innovation in molecular imaging and radiopharmaceutical therapy. Previously, he served as senior attending physician at the Theranostics Center in Bad Berka, Germany—one of the world's pioneering sites for PSMA-, SSTR-, and FAP-targeted care—where he was instrumental in advancing global practice in theranostics and molecular imaging. Respected for his thought leadership and scientific insight, Dr. Kulkarni has authored numerous publications, presented widely at international fora, and remains committed to expanding access to advanced theranostics worldwide.

PASSION for CHANGE

©2025 United Imaging Healthcare Co., Ltd. All rights reserved.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you are welcome to contact us at the following email address: compliance@united-imaging.com