

u!nnovation.

**UNITED
IMAGING** 



Global Spread of external contributors



Foreword

A Note from the Editorial Board

Welcome to this year's annual magazine from the Global Scientific Collaboration team.

uInnovation is a scientific magazine published by United Imaging Healthcare that has been successfully distributed for over past three years. It aims to serve as a platform for sharing ground-breaking advancements, emerging trends, and future possibilities in the vast expanse that is oncology.

uInnovation is currently in its fourth edition. This year's edition will inform, engage, and inspire you about the latest developments and applications of United Imaging Healthcare. This journal includes quick read sections for those in a rush, and appealing images to promote visual understanding.

United Imaging Healthcare is a provider of high-end medical equipment and medical IT solutions. From our headquarters in Shanghai's Jiading district to our network of research and development centers throughout the world, our global mission is to provide medical institutions with a full range of healthcare solutions, from diagnostic imaging and radiation therapy equipment to service, training, and medical IT solutions.

We are dedicated to expanding access to high-quality medical care and improving the value of our services.

VISION — Leading Healthcare Innovation

MISSION — To Bring Equal Healthcare for All

To learn more, visit <https://www.united-imaging.com/>

Disclaimer

The articles contained in this magazine are provided solely by the authors, and the author(s) of each article appearing in this magazine is/are solely responsible for the content thereof as well as personal data, which is used anonymously or complied with applicable data privacy laws or regulations. United Imaging Healthcare makes no representation or warranties, expressly or impliedly, with respect to the accuracy, timeliness, reliability, legitimacy, applicability, fitness, originality, or completeness of the contents of this magazine. United Imaging Healthcare assumes no legal responsibility or liability for any error, omission, or illegality with respect to the material contained within.

All articles contained in this magazine only represent the opinions and views of the authors and do not implicitly or explicitly represent any official positions or policies, or medical opinions of United Imaging Healthcare or the institutions with which the authors are affiliated unless this is clearly specified. Discussions of any brand, services, or products in the magazine should not be construed as promotion or endorsement thereof.

Articles published in this magazine are intended to inspire further general scientific research, investigation, understanding, and discussion only and are NOT intended to and should not be relied upon as recommending or promoting a specific medical advice, method, diagnosis, or treatment by physicians for any particular individual, nor to replace the advice of a medical doctor or other healthcare professional. Any individual wishing to apply the information in this magazine for the purposes of improving their own health should not do so without consulting with a qualified medical practitioner. All patients need to be treated in an individual manner by their personal medical advisors. The decision to utilize any information in this magazine is ultimately at the sole discretion of the reader, who assumes full responsibility for any and all consequences arising from such a decision. United Imaging Healthcare makes no representations or warranties with respect to any treatment, action, or application of medication or preparation by any person following the information offered or provided within or through the magazine. United Imaging Healthcare shall remain free of any fault, liability, or responsibility for any loss or harm, whether real or perceived, resulting from the use of information in this magazine.

The articles included in this magazine may contain work in progress, which represents ongoing research and development. Such technologies are not available for sale in the United States for clinical use and also may not be available for such sales in other countries around the world.

Please note that the magazine is intended to be distributed only within a limited scope instead of publication.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you are welcomed to contact us at the following email address: compliance@united-imaging.com

Clinical application of CT-Linac with Artificial Intelligence for Online Adaptive Radiotherapy in Head & Neck and Pelvic Cancers

Guanqun Zhou^{*a}, Shuai Sun^{*b}, Lecheng Jia

^a Department of Radiation Oncology, Sun Yat-sen University Cancer Center, State Key Laboratory of Oncology in South China, Collaborative Innovation Center for Cancer Medicine, Guangdong Key Laboratory of Nasopharyngeal Carcinoma Diagnosis and Therapy, Guangzhou, CHINA.

^b Department of Radiation Oncology, Peking Union Medical College Hospital, Chinese Academy of Medical Sciences & Peking Union Medical College, Beijing, CHINA.

*Shared Authorship

The Challenge

"Clinics struggle with inefficient workflows. Traditional offline replanning takes days. What's missing is a fast, accurate, daily adaptive solution that fits into busy clinics without overwhelming staff or extending treatment time"

Delivering precise radiation to tumors is complicated by the fact that a patient's anatomy changes daily — organs shift, tumors shrink, and weight fluctuates. This is especially problematic in cancers like nasopharyngeal and cervical cancer, where targets sit close to sensitive organs. When the original treatment plan no longer matches the patient's current anatomy, tumors may be underdosed, or healthy tissues overdosed leading to poor control or unnecessary side effects like hearing loss, xerostomia, or bowel toxicity.

The Solution

The uRT-Linac 506c — a diagnostic-quality fan-beam CT (FBCT)

scanner fully integrated with a linear accelerator, developed by United Imaging Healthcare (UIH) offering these features:

Precision Imaging by FBCT: Same quality as CT Sim, can be used for contouring and dose calculation directly.

All-in-One Platform: Simulation, daily imaging, contouring, planning, and delivery happen on one machine — no patient transfer, no image mismatch.

AI Auto-Segmentation: Uses deep learning to contour targets and 35+ organs in nearly one minute, editable by clinicians.

AI Auto-Planning: Generates deliverable VMAT plans using prior plan knowledge and clinical goals less than 3.5 minutes (nasopharyngeal cancer) and less than 2.5 minutes (cervical cancer).

Real-Time QA: In-vivo EPID dose monitoring with gamma analysis during beam-on.



Figure 1. The integrated uRT-Linac 506c from United Imaging Healthcare #.

[#]No 510k application for uRT-Linac 506c has been filed with the FDA. This product is not available for sale in the U.S. for clinical uses and also may not be available for sales in other countries.





Figure 2. The comparison of simulation CT and FBCT by CT-Linac of a representative patient. The top row is simulation CT and the bottom row is FBCT.

Workflow Summary

The online adaptive radiotherapy (oART) workflow using CT-Linac consists of the following steps:



Did that Work?

FBCT-guided oART demonstrates clear dosimetric advantages in definitive cervical cancer treatment by ensuring consistent, high-quality target coverage while

This turns oART from a research-only setting into a clinically feasible daily routine.

Here are the average RT times for two different types of cancer:

- i. For Nasopharyngeal carcinoma, the following are the average time fraction:
 - FBCT Scan & Registration: ~3.3 min (16% of time)
 - AI Contouring + Review: ~8.2 min (39% of time)
 - AI Planning + Approve: ~5.0 min (24% of time)
 - Beam Delivery + EPID in Vivo: ~4.4 min (21% of time)
- ii. For cervical cancer, the following are the average time fraction:
 - FBCT Scan & Registration: ~2.3 min (10% of time)
 - AI Contouring + Review: ~10.1 min (44% of time)
 - AI Planning + Approve: ~5.0 min (22% of time)
 - Low-dose verification FBCT: ~3 min (13% of time)
 - Beam Delivery + EPID in Vivo: ~2.4 min (11% of time)

This streamlined workflow demonstrates that CT-Linac-based oART can be efficiently integrated into routine clinical practice with acceptable time and dosimetry accuracy.

significantly reducing radiation exposure to critical organs at risk—thereby translating anatomical adaptability into tangible improvements in plan accuracy and safety.

Representative data from Nasopharyngeal carcinoma:




Figure 3. Dose distribution comparison between (a-c) scheduled plans and (d-f) ART plans. (adopted from [1])




Figure 4. Dose and Volume Histogram (DVH) of targets(a) and OARs(b) of scheduled plans (dotted line) and ART plans (solid line). (adopted from [1])

The cumulative DVH plot of targets and OARs are shown in Figure 4, respectively. The solid line represents the ART plan, and the dotted line represents the scheduled plan.

ART plan is better than the scheduled plan - the dose to critical organs is reduced in ART. For example, in chiasm, there is 43% lower dose in ART plan.

Representative data from a Cervical cancer case:

Dose volume histogram comparison are shown for a cervical cancer patient. Adaptive planning consistently improved target coverage (PTV D95 met prescription) while reducing

dose to rectum, bladder, and small bowel in Figure 5. Real-time EPID monitoring confirmed >99% gamma passing rate, validating delivery accuracy. Even with tumor movement and bladder fullness, adaptive plans-maintained precision — something static plans failed to do.




Figure 5. The differences in dosimetric parameters of targets and OARs between ART plans(red) and scheduled plans(blue) : (A) Targets (B) OARs. (adopted from [2])

The solid line represents the ART plan, and the dotted line represents the scheduled plan. So, it shows that ART plan is better than Scheduled plan because of improved distribution to both targets and organs at risk (OARs). First, the dose is reduced in ART plans for critical OARs. For example, in the small intestine, there is a 2.17 Gy lower D2cc (dose to the most irradiated 2 cm³) for the ART plan, and in the rectum, the D2cc is 0.10 Gy lower, demonstrating superior organ sparing.

So What?

Patients win: fewer side effects, better tumor targeting, preserved organ function (like ovaries or hearing). Clinics win: no workflow disruption, no extra machines, no long waits.

The future? Daily adaptation becomes standard — not for select cases, but for all patients who need precision.

UIH's platform sets a new benchmark: adaptive radiotherapy that's fast, smart, and ready for prime time.

The impact? Higher cure rates, better quality of life, and a scalable model for global implementation — even in resource-limited settings.




Figure 6. The DVH of the scheduled plans(dotted line) and ART plans(solid line) for a representative patient (adopted from [2]) .


References:

- [1] Yang, Y-X, Li L et al. (2024). Initial Experience of CT-Based Online Adaptive Radiotherapy for Nasopharyngeal Carcinoma With a Novel Integrated Platform: A Case Report. Cureus.
- [2] Sun, Gong et al (2025). Evaluating the Implementation of fan-beam CT-guided online adaptive re-planning in definitive cervical cancer radiotherapy. Frontiers in Oncology.

Author Statement

Until now, adaptive radiotherapy was too slow, too complex, or too inaccurate for routine use. Our integrated CT-Linac with AI automation makes daily adaptation not just possible — but practical for busy departments treating complex head and neck or pelvic cancers.

Author’s Biography



Dr. Guanqun Zhou
Department of Radiation Oncology, Sun Yat-sen University Cancer
Centre, State Key Laboratory of Oncology in South China,
Collaborative Innovation Centre for Cancer Medicine, Guangdong
Key Laboratory of Nasopharyngeal Carcinoma Diagnosis and
Therapy, Guangzhou, CHINA

Dr. Guanqun Zhou, MD, PhD, is a radiation oncologist at the Department of Radiation Oncology, Sun Yat-sen University Cancer Center (SYSUCC) in Guangzhou, China. She specializes in treating head and neck cancers, particularly nasopharyngeal carcinoma. Dr. Zhou is also active in academic research, leading work on online adaptive radiotherapy (ART) using integrated CT-guided platforms, showing feasibility and dosimetric benefits for NPC. Her research contributions include studies on AI-driven planning, deep learning segmentation of tumor volumes, and the dynamics of tumor/organ changes during radiotherapy.



Dr. Shuai Sun
Department of Radiation Oncology,
Peking Union Medical College
Hospital, Chinese Academy of Medical
Sciences & Peking Union Medical
College, Beijing, CHINA

Dr. Shuai Sun is a radiation oncologist at the Department of Radiation Oncology, Peking Union Medical College Hospital, where he is involved in clinical care, research, and the advancement of precision radiotherapy. Her academic work focuses on modern radiation therapy techniques—including IMRT, VMAT, stereotactic radiotherapy, and image-guided treatments—with particular interest in applying quantitative imaging, treatment-response assessment, and technology-driven workflow optimization in oncology. Dr. Sun has contributed to peer-reviewed publications and collaborative clinical studies addressing treatment outcomes and toxicity reduction, and he actively participates in multidisciplinary cancer management at one of China’s leading national medical centers.

PASSION for CHANGE

©2025 United Imaging Healthcare Co., Ltd. All rights reserved.

If you have any questions about the magazine, or simply wish to reach out to us for any other reasons, you are welcome to contact us at the following email address: compliance@united-imaging.com